CD34Count Kit

Code No. K2370/ Réf. K2370/ Code Nr. K2370

1st edition/ 1ère édition/ 1. Ausgabe

For the identification and enumeration of CD34+ cells in human mobilized blood samples and leukapheresis samples.

The kit contains reagents sufficient for 50 duplicate tests.

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Page/ Page/ Seite

ENGLISH

Intended Use	5
Summary and Explanation	5
Principle of the Procedure.....	5
Reagents	6
Materials provided	6
Materials required but not provided	6
Reagent Preparation.....	7
Precautions.....	7
Storage	8
Specimen Collection and Preparation	8
Staining Procedure	8
Staining of cells and addition of count control beads	8
Guideline, “wet tip” reverse pipetting technique.....	9
Quality Control	10
Control sample	10
Time as quality control	10
Duplicate tests.....	10
Data Acquisition	10
Acquisition set-up for CD34Count Kit	10
Analysis of Data	13
Results.....	13
Value Assignment.....	14
Limitations	14
Troubleshooting	15
Instrument.....	15
References	40
Explanation of symbols	40

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ENGLISH

Intended Use

For in vitro diagnostic use in flow cytometry. CD34Count Kit is intended for use as an in vitro diagnostic test to identify and enumerate CD34-positive (CD34+) haematopoietic stem cells in human mobilized peripheral blood samples and leukapheresis samples.

Summary and Explanation

Haematopoietic progenitor cells (HPCs) can be mobilized from the bone marrow into the peripheral blood (PB) by cytotoxic drugs, cytokines, or combinations of the two. This mobilization allows the collection of HPCs by apheresis in sufficient quantities for transplantation procedures (1, 2). The use of PB stem cells (PBSCs) for transplant purposes has rapidly increased during the past decade, and for autologous transplantation PB is now preferred as a source of HPCs over bone marrow (BM). Traditionally, the absolute mononuclear cell count of BM harvests in relation to the body weight of the recipient has served as a useful predictor of engraftment potential. However, this approach is insufficient for mobilized PBSC collections owing to their variable HPC content. By the late 1980s, it was established that virtually all of the CFC (colony-forming cell) activity and engraftment potential of marrow or peripheral blood samples were contained in the small population of cells bearing the CD34 antigen (3). CD34 is a cell surface antigen, which on haematopoietic cells is restricted to early progenitors of all lineages (1, 2). Thus, CD34+ HPCs can restore multilineage haematopoiesis in myelo-ablated patients. The number of CD34+ cells as determined by flow cytometry per kilogram of recipient body weight has been shown to be the most useful indicator of the haematopoietic reconstitutive capacity of peripheral blood stem cell (PBSC) transplants (1, 2).

Sutherland et al. (5) developed a set of clinical guidelines for the "International Society for Hematotherapy and Graft Engineering" (ISHAGE) for simple and standardized CD34+ enumeration. These guidelines have been adopted in the CD34Count Kit that contains Anti-CD45/FITC and Anti-CD34/RPE, and thus allows the identification of leukocytes (CD45+) and the verification of "true" CD34+ cells as being dim for CD45 fluorescence and having low side scatter (CD45^{dim}, SSC^{low}). The addition of the CytoCount™ counting beads to the cell suspension yields the concentration of CD34+ cells per unit of original specimen volume (i.e. the absolute CD34+ cell count) from a single flow-cytometric assessment (single-platform technique) (1-3). The single-platform technique has fewer sources of variability compared with the previously commonly used dual-platform technique, and multicentre studies have shown that it provides a better coefficient of variation (CV), and that the risk of generating aberrant results is lower (4).

The optional addition of the DNA stain, 7-aminoactinomycin D (7-AAD), which is also included in the kit, makes it possible to exclude dead cells from analysis (1, 2). 7-AAD is only able to diffuse into and stain dead cells, because the cell membrane in these cells is no longer intact.

Principle of the Procedure

CD34Count Kit has been designed to provide an optimal method for CD34+ cell enumeration, following the ISHAGE guidelines for single-platform (1, 2, 5).

The test is performed in duplicate by staining the sample with a dual-colour reagent composed of Anti-Human CD45/FITC and Anti-Human CD34/RPE. After staining of the sample, red blood cells are lysed with an ammonium chloride-based lysing reagent, and 7-AAD is added (optional) to the sample if viability measurement is wanted. Fluorescent count control beads (CytoCount™) are added to the sample for use as an internal reference population.

During analysis, the absolute concentration of CD34+ cells in the sample can be determined by dividing the number of CD34 events by the number of fluorescent bead events, and then multiplying by the bead concentration in CytoCount™.

The principle of the ISHAGE guidelines is to use sequential Boolean gating strategy to define the real progenitor cells with regenerative potential. This is obtained by counter-staining CD34+ cells with Anti-CD45/FITC and thus allowing the elimination of debris and non-specifically stained events from the analysis (1, 2). Importantly, this approach also allows the discrimination between HPCs, which express relatively low levels of CD45 on their surface, and lymphocytes and monocytes, which express high levels. Just as lymphocytes, monocytes, and granulocytes form discrete clusters on bivariate plots of CD45 versus side scatter (SSC), so do non-malignant CD34+ haematopoietic stem and precursor cells (1, 2). The addition of a known concentration of reference beads makes it possible to calculate the absolute concentration of CD34+ cells in the original specimen, i.e. to determine the number of CD34+ cells per μL of specimen.

Exclusion of dead cells from CD34+ cell enumeration is optional and can be achieved by inclusion of the DNA intercalating dye, 7-AAD, which allows discrimination between viable and nonviable cells. 7-AAD is excited at 488 nm and has emission maximum at 660 nm. Addition of 7-AAD is particularly recommended when analysing leukapheresis samples stored for more than four hours or leukapheresis samples analysed after thawing.

Reagents

Materials provided

Vial 1

ANTI-HUMAN CD45/FITC
ANTI-HUMAN CD34/RPE

**Monoclonal Mouse Anti-Human CD45/FITC, Clone T29/33 +
Monoclonal Mouse Anti-Human CD34/RPE, Clone BIRMA-K3**
1 mL, ready to use

Vial 2

EasyLyse™
2 x 5 mL, 20 x concentrated
Erythrocyte-lysing reagent, ammonium chloride based

Vial 3

7-AAD

7-Aminoactinomycin D (7-AAD)
1 mL (0.01% w/v), ready-to-use
For viability staining

Vial 4

CytoCount™
17 mL, ready-to-use after resuspension
Count control beads

Materials required but not provided

1. 12 x 75 mm polystyrene tubes.
2. Calibrated pipette for 10 μL and 100 μL volume.
3. Blood mixer or rotary mixer, for re-suspension of the CytoCount™ beads.
4. Deionized or distilled water.
5. Vortex mixer.
6. Bulk dispenser or pipettor (2 mL) for dispensing 2 mL EasyLyse™.
7. PBS for dilution of specimen if necessary.
8. Sheat Fluid (DakoCytomation code No. S2322).
9. CD34+ control samples (DakoCytomation CD-Chex® CD34, code No. K2358 or K2422).

Reagent Preparation

1. Dilute the required amount of EasyLyse™, Vial 2, 1+19 with deionized water of room temperature. Alternatively mix the contents of one bottle (5 mL) of EasyLyse™ with 95 mL of deionized water of room temperature. Tap water is not compatible with the lysing reagent.
2. Remove Vial 4, CytoCount™, from the 2-8 °C storage. Make sure the cap is tightened to prevent fluid leakage and mix well to ensure complete re-suspension of the beads. Avoid mixing with a vortex mixer owing to the potential introduction of air bubbles that may reduce the actual volume of beads pipetted. It is advisable to place the beads on a rotary mixer for 5-60 minutes to ensure that the beads remain in suspension until required. Allow the CytoCount™ reagent to reach room temperature before use.

Precautions

1. For professional users.
2. Vial 1, Anti-Human CD45/FITC + Anti-Human CD34/RPE, contains material of animal origin. As with any product derived from biological sources, proper handling procedures should be used.
3. Do not freeze any of the reagents or expose any of the reagents to direct light.
4. Vial 3, 7-AAD, is a nucleic acid-binding dye. The concentration of dye in the solution is 0.01% and therefore the solution does not require hazard labelling. The toxicological properties of this dye have not been investigated.
5. The concentration of beads in Vial 4, CytoCount™, varies from lot to lot. The bead concentration (C) and standard deviation (S) are stated on the vial as number of beads per μL .
6. Vials 1, 2, 3 and 4 contain sodium azide (NaN_3), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
7. Ensure that Vial 4, CytoCount™, is stored tightly capped in the upright position to prevent fluid leakage. The CytoCount™ beads must be thoroughly re-suspended before use. Harsh mixing of the CytoCount™ beads on a vortex mixer can introduce air bubbles into the fluid, and this can reduce the expected volume of CytoCount™ beads pipetted. Gentle rotation (e.g. on a rotary mixer or blood tube roller) is recommended.
8. Care must be taken not to lose any of the CytoCount™ suspension before re-suspension; otherwise the bead count becomes invalid. If any leakage from Vial 4, CytoCount™, has occurred upon receiving the product, do not use the beads. Contact DakoCytomations Technical Service for replacement of Vial 4.
9. It is essential to use the same calibrated pipette and the "wet tip" reverse pipetting technique for both the sample and the CytoCount™ beads. It is important that the pipetted volume of the CytoCount™ beads and the sample is accurate. The volumes of the antibody reagent and lysing reagent are not critical in connection with determination of the absolute CD34+ cell count number; it is the ratio between the volume of the sample and the volume of the beads that is important.
10. It is recommended to set the threshold value on FL1 (515-545 nm), since this minimizes the risk of losing CytoCount™ events during acquisition. Care should be taken not to exclude any CD34+ events by setting the FL1 threshold too high. Please refer to Troubleshooting.
11. Testing showed no carry over between two different samples. It is however recommended that each flow cytometer is validated for carry over effect between samples. To avoid any carry over from a sample with a high to a sample with a lower CD34+ cell concentration, PBS can be run through the flow cytometer between samples.
12. Optimal compensation settings are critical when using 7-AAD as a marker for dead cells because the signal from the RPE-channel overlaps into the 7-AAD-channel. In the worst-

case scenario, positive CD34 events will be excluded if they fall into the dead-cell exclusion gate in the 7-AAD vs. Side Scatter (SSC) channel.

13. 7-AAD is excited at 488 nm and has emission maximum at 660 nm (FL3 or FL4, depending on the instrument). Consequently, the dye cannot be used in single-laser systems with other fluorochromes that emit at >600 nm, such as PerCP or RPE-Cy5.
14. Aggregation of count control beads has been reported (4). CytoCount™ beads exhibit an especially low number of multiples (below 3%) making the gating of the beads more simple. If a higher number of CytoCount™ multiples is observed and a problem with the product is suspected, please contact DakoCytomations Technical Services.

Storage

Store the kit at 2-8 °C in the dark. Do not use after expiration date stamped on the kit. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user.

EasyLyse™, Vial 2, is stable for 6 months after opening of the vial. After dilution the lysing reagent is stable for 20 days at 2-24 °C.

Store all reagents in the upright position.

Specimen Collection and Preparation

Collect peripheral blood aseptically by venipuncture into a sterile K₃EDTA blood collection tube. Each test requires 200 µL of sample (100 µL per tube). Follow the collection tube manufacturer's guidelines for sample collection. Collect leukapheresis samples according to the instructions from the manufacturer of the leukapheresis instrument.

Store anti-coagulated blood at room temperature (20 to 25 °C) until ready for staining.

Stain peripheral blood samples within 24 hours of collection. It is recommended to stain leukapheresis samples at the time the leukapheresis product is being processed.

Perform a white blood cell (WBC) count on all specimens to be evaluated. If the WBC count is greater than 10⁷ WBC/mL, dilute the specimen to about 10⁷ WBC/mL with PBS. Record the Dilution Factor for the calculation of the final absolute CD34+ cell count in the original specimen.

Do not use pre-fixed cells.

Staining Procedure

Staining of cells and addition of count control beads

1. Dilute the original specimen to about 10⁷ leucocytes per mL with PBS. Record the accurate Dilution Factor. In most cases, dilution will be required only for leukapheresis samples.
2. Every sample must be run in duplicate. Therefore mark two test tubes for every sample.
3. Pipette 100 µL of the sample into each test tube. Use reverse pipetting (see Guideline, "wet tip" reverse pipetting technique below).
4. Add 10 µL of Vial 1 (Anti-Human CD45/FITC + Anti-Human CD34/RPE) to each tube and mix on a vortex mixer. Incubate the tubes for 30 minutes at 2-8 °C or for 15-30 minutes at room temperature (20-25 °C).
5. Add 2 mL of diluted lysing reagent (Vial 2, EasyLyse™ diluted 20 x). Mix the sample on a vortex mixer immediately after adding the lysing reagent.
6. Incubate in the dark for 10 minutes at room temperature.
7. If viability staining is wanted add 10 µL of Vial 3 (7-AAD) to each sample and mix. Otherwise proceed to step 8. Incubate at least 5 min with 7-AAD before acquisition takes place.
8. Resuspend Vial 4 (CytoCount™) on a rotary mixer as described in Reagent Preparation. Add 100 µL to the stained and lysed sample. Pipetting of both the sample and CytoCount™

beads should be performed with the same calibrated pipette and using a “wet tip” reverse pipetting technique (4) (see Guideline, “wet tip” reverse pipetting technique below).

9. Mix beads and sample gently using a vortex mixer for 3 seconds before acquisition to ensure an even distribution of CytoCount™ beads throughout the sample.
10. Sample acquisition must take place within 45 minutes after the addition of the lysing reagent (step 5).
11. Acquire a minimum of 60 000 white blood cell (WBC) events, 1000 CytoCount™ beads, and a minimum of 100 CD34+ cells (4). Please refer to Acquisition set-up for the CD34Count Kit.

Guideline, “wet tip” reverse pipetting technique

The same pipette should be used to dispense sample and beads. Use a clean pipette tip for each reagent. Adopting the same technique for dispensing both sample and beads will ensure more consistent results.

Sample

1. Press plunger to the second stop, place tip in the fluid and aspirate the sample. Then press the plunger to the first stop dispensing the sample. The pipette tip should now contain excess fluid.
2. Repeat this step twice keeping the pipette tip within the sample.
3. The pipette tip is now ready for use. Press plunger to the **second** stop and place tip in the fluid and aspirate the sample. Carefully remove the pipette tip from the sample and wipe the pipette tip gently with tissue to remove any excess fluid, being careful not to touch the tip orifice and inadvertently removing sample from inside the pipette tip.
4. Try to keep the pipette as vertical as possible. Check for the presence of air bubbles in the tip. If air bubbles are present then dispense the fluid back into the sample and repeat the sampling procedure.
5. Dispense the sample to the **first** plunger stop, residual fluid should be seen in the pipette tip after the fluid has been dispensed. Do not touch the walls of the test tube when removing the tip.

Beads

6. Prepare the pipette with the CytoCount™ bead suspension in Vial 4 as described above for the sample (steps 1-3).
7. For dispensing bead suspension, place the pipette tip on the tube wall near to the top of the sample; slightly angling the tube and pipette will aid this step. Dispense to the **first** plunger stop, residual fluid should be seen in the pipette tip after the fluid has been dispensed.
8. For further pipetting of beads, place the tip in the CytoCount™ Vial and aspirate more bead suspension without dispensing the residual fluid left in the tip after the first pipetting.
9. If the tip has not been in touch with the sample, the residual fluid in the tip should be dispensed into the CytoCount™ Vial. If you choose to discard the residual fluid, there might no longer be reagent for 50 duplicate tests.

Quality Control

Control sample

It is recommended to include a positive control sample as quality control check for a proper CD34+ cell enumeration set-up. DakoCytomation CD-Chex[®] CD34, code No. K2358 or K2422 is recommended as positive control sample. When using CD-Chex[®] CD34, 7-AAD should not be added to the sample, since 7-AAD will stain all cells because the control cells have been fixed. Participation in an external quality assessment scheme is highly recommended.

Time as quality control

Time can be used as an internal quality control check for single-platform determinations, since studies have demonstrated that the number of beads acquired in a fixed time period is remarkably constant for most flow cytometers. Thus for each new batch of beads, the user can derive a mean bead count ± 2 SD range for a given time. Any sample with a bead count outside of this range is indicative of a potential pipetting error and the sample should be re-stained (6).

Duplicate tests

Each specimen must be tested in duplicate. The difference between the values obtained by the duplicate tests should not exceed 15%. Otherwise the sample should be re-tested. The difference between the duplicate values (X and Y) is calculated as follows:

$$\frac{|X-Y|}{(X+Y)/2} \times 100\%$$

Data Acquisition

Follow the manufacturer's recommendations for setting up the cytometer for three-colour acquisition. Set-up should include setting or checking photomultiplier tube (PMT) voltages and instrument sensitivity, as well as compensation, if hardware compensation is used.

Acquisition set-up on FACSCalibur (Beckton Dickinson) and on the Coulter XL (Beckman Coulter) flow cytometers has been described by Sutherland et al. (1) and Gratama et al. (2).

Acquisition set-up for CD34Count Kit

The data can either be acquired with hardware compensation while acquisition takes place, **or** they can be acquired without compensation, which will then be applied during analysis of the data. Please refer to common flow cytometry guidelines regarding compensation set-up.

The following protocol is designed both for acquisition and for analysis.

1. Generate the following two-dimensional data displays (dot plots) and set regions according to the figures:

Plot A: CD45 FITC vs. SSC-H (Region 1 and Region 5)

Plot B: CD34 RPE vs. SSC-H (Region 2)

Plot C: CD45 FITC vs. SSC-H (Region 3)

Plot D: FSC-H vs. SSC-H (Region 4)

Plot E: 7-AAD vs. SSC-H (Region 8)

Plot F: CD45 FITC vs CD34 RPE (Region 6 and Region 9)

Plot G: Time vs. FSC-H gated on the CytoCount™ events (Region 7)

Plot H: FSC-H vs. SSC-H

2. Create the following gate logic (colour is optional):

Name	Expression	Color
G1	NOT R8	
G2	R1 & NOT R8	Green 16 Shades
G3	R1 & R2 & NOT R8	
G4	R1 & R2 & R3 & NOT R8	
G5	R5 & NOT R8	
G6	R6	
Beads	R6 & R7	Blue 16 Shades
Stem Cells	R1 & R2 & R3 & R4 & NOT R8	Red 16 Shades

G7 is named “Stem Cells” and is added in the gate logic builder to colour the positive CD34+ events. The colour gating is very useful during acquisition set-up and analysis.

3. Apply the gates onto the plots as follows

Plot A: G1

Plot B: G2

Plot C: G3

Plot D: G4

Plot E: Ungated

Plot F: Ungated


Plot G: G6

Plot H: Ungated

4. Set the threshold value on FL-1 so that it excludes all CD45-negative events. It is recommended to set the threshold value on FL-1, since this minimizes the risk of exclusion of the CytoCount™ beads. Threshold setting on the FSC parameter is feasible, but care should be taken not to exclude any bead events. Please refer to Troubleshooting

5. Load the first sample and adjust the PMT gain for FSC and SSC by looking at Plot H to ensure all populations are on scale. If the threshold is set on the FSC parameter, make sure that the CytoCount™ events are placed over the threshold value.


6. Adjust the PMT for CD45/FITC. The lymphocyte population should fall within the third decade (Plot A). Make sure that no dim CD45+ cells are excluded. If in doubt, either lower the FL1 threshold or raise the PMT gain for FL1. If PMT settings are changed, it might be necessary to adjust the compensation settings as well.

7. Adjust the PMT for CD34/RPE. The majority of the negative events should fall between the first and second decade (Plot B).


8. Adjust the PMT for 7-AAD. The majority of the negative events should fall within the first decade (Plot E).

9. Apply Region 6 from Plot F in Plot G. If the threshold is set on the FSC parameter, use Plot G to control that no CytoCount™ events have been excluded owing to the threshold setting (refer to Troubleshooting).

10. It is recommended to apply an exclusion gate in Plot H. This exclusion gate shall include the SSC-low and FSC-low events (Region 9). During acquisition the events falling into this gate should be rejected. This is obtained by applying a gate filter (exclusion gate) for these events or by rejecting events in R9 in the acquisition window. When the FSC parameter is used as threshold this exclusion gate is not relevant.

11. Make sure that time is saved as a parameter during acquisition.

The instrument is now ready for acquisition.

12. Acquire at least 60 000 WBC events, 1000 CytoCount™ events and 100 CD34+ events. It is only possible to set a stop criterion for one of these numbers. For low CD34+ samples the 100 CD34+ events will be the stop criterion, while the 60 000 WBC events will be the stop criterion for high CD34+ samples. It is advisable to acquire as many events as feasible, since this will provide the most accurate CD34+ cell number.

Analysis of Data

Use the same protocol for analysis as used for acquisition.

1. Apply compensation setting by following the software instructions for data acquired without hardware compensation. Note that compensation for overlap from the RPE channel into the 7-AAD channel is critical (refer to Precautions).
2. Statistics shall be shown for Region 4 (Stem Cells) showing CD34+ progenitor cells, and for Region 7 (Beads) showing the bead count. Statistics should also be shown for Region 1 (G2), if the number for acquired living CD45-positive cells is required.
3. Import the acquired data file (FCS file) and position the regions.
4. In Plot D apply G5 and position Region 4 to include only the smallest lymphocytes (Figure 1(a)).
5. Apply now G4 in Plot D and keep the Region 4 in the same position (Figure 1(b)).


Figure 1. (a): Region 4 is positioned to include events bigger than the smallest lymphocytes as well as events low in Side Scatter. (b): after setting Region 4 on living lymphocytes, the plot is gated on G4 to show those events fulfilling the requirements described in the ISHAGE gating strategy.

Results

Calculation of the absolute CD34+ cell count per μL of the original specimen:

Use the following equation to calculate the absolute CD34+ cell count:

CD34+ Cells/ μL =

$$\frac{\text{Counted Number of CD34+ Cells (R4)} \times \text{Concentration of CytoCount}^{\text{TM}*}}{\text{Counted Number of Beads (R7)}} \times \text{Dilution Factor}$$

*This value (C) is found on the label of Vial 4: C = Number of beads/ μL .

In the example shown (Plots A to H), the absolute CD34+ cell count calculation is as follows:

$$\frac{216 \times 1006}{3276} \times 1 = 66 \text{ CD34+ Cells}/\mu\text{L}$$

The final result is the average of the duplicate tests. DakoCytomation recommends that the difference between these two values does not exceed 15%, otherwise the sample should be analysed again (refer to Quality control).





Value Assignment

The bead concentration (C) and corresponding standard deviation (S) is stated on the label of Vial 4 as number/ μL . This number is obtained by using a Coulter Particle Counter Z1. The bead concentration standard deviation is a minor contributor to the combined uncertainty of the analytical result.

Limitations

1. Because of the variability of biological samples, some samples might not be correctly analyzed. An experienced flow cytometrist must critically review all dot plots. An incorrectly high result could lead to an infusate with less than the recommended threshold dose of CD34+ cells.
2. CD34Count Kit results must be interpreted in conjunction with other information available from clinical evaluation and additional diagnostic procedures.
3. Store blood samples at room temperature (20 to 25 °C).
4. Do not use previously fixed and stored patient specimens.
5. The WBC count should not exceed 10^7 WBC/mL and the CD34 cell count should not exceed 2000 cells/ μL in the sample.

Troubleshooting

Problem	Solution
<p>Some bead events are lost because the threshold value on the FSC-H parameter has been set too high.</p>  <p>Wrong: In the FSC-H parameter no events are observed beneath the main bead population. Some CytoCount™ beads have been excluded during acquisition.</p>	<p>Re-run the sample with a lower threshold on FSC-H. However it is recommended to set the threshold on the FL1 parameter.</p>  <p>Correct: Events smaller than the main bead population can be seen. CytoCount™ beads have not been excluded.</p>
<p>Part of the CD34+ population is excluded because the FL1 threshold has been set too high.</p> 	<p>Re-run the sample with a lower FL1 threshold.</p> 
<p>Quality Control Sample is out of range.</p>	<p>Check your flow cytometer and acquisition set-up, and run the Quality Control Sample again. If the control sample is still out of range, contact your local instrument support.</p>
<p>It is difficult to set the regions.</p>	<p>Enlarge the dot plots and reposition the region.</p>

Instrument

CD34Count Kit is intended for use on a flow cytometer. The flow cytometer must be equipped with appropriate computer hardware and software and with a 488 nm laser. The flow cytometer must be capable of detecting fluorescence at the following wavelengths: 515-545 nm (FL1), 562-607 nm (FL2) and >650 nm (FL3 or FL4, depending on the instrument). Also, light scatter (forward and side scatter) must be detectable. The instrument must be able to threshold on either the FL1 or the forward scatter parameter.

Follow the quality control set-up recommended by the instrument manufacturer. Set-up should include setting or checking PMT voltages, fluorescence compensation, and instrument sensitivity.

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1. Sutherland DR, Keeney M, Gratama JW. Enumeration of CD34+ hematopoietic stem and progenitor cells. Current Protocols in Cytometry. New York: John Wiley & Sons; 2003. p. 6.4.1-23.
2. Gratama JW, Keeney M, Sutherland DR. Enumeration of CD34+ hematopoietic stem and progenitor cells. Current protocols in cytometry. New York: John Wiley & Sons; 1999. p. 6.4.1-22.
3. Keeney M, Gratama JW, Sutherland DR. Critical role of flow cytometry in evaluating peripheral blood hematopoietic stem cell grafts. Cytometry 2004;58A:72-5.
4. Brando B, Barnett D, Janossy G, Mandy F, Autran B, Rothe G, et al. Cytofluorometric methods for assessing absolute numbers of cell subsets in blood (review). Cytometry 2000;42:327-46.
5. Sutherland DR, Anderson L, Keeney M, Nayar R, Chin-Yee I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematother 1996;5:213-26.
6. Bergeron M, Lustyik G, Phaneuf S, Ding T, Nicholson JKA, Janossy G, et al. Stability of currently used cytometers facilitates the identification of pipetting errors and their volumetric operation: "time" can tell all. Cytometry 2003;52B:37-40.

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DakoCytomation Denmark A/S